



EIE INSTRUMENTS PVT. LTD.
(An ISO 9001:2015 Certified Company)

SAND BATH-DIGITAL-30X20X10 CM



Applications :

A laboratory sand bath is a, generally, electrically heated container filled with sand used for uniform, gradual heating of glassware, such as flasks or beakers, without direct heat contact. Offering superior temperature stability, they are ideal for drying, distilling, and heating sensitive samples upto 300 °C

Salient Features :

1. Aesthetic outer appearance
2. Heated by open type heaters placed at the bottom of the inner chamber
3. 3 stage heat control switch for selection of LOW, MEDIUM or HIGH heat
4. Fully air ventilated lower body for minimum heat, even when the top plate is fully heated
5. Uniform heating throughout the bottom Plate

Optional Accessories At extra cost :

1. Forced air exhaust for continuous operation
2. Thermal gloves for thermal protection.

Notes:

- Thermal Inertia: Sand baths stay hot for a significant amount of time after the heater is turned off.
- Slow Warm-up: Sand has low thermal conductivity, so it can take 90 minutes or longer to reach a stable high temperature.
- High Thermal Inertia: Once heated to its maximum capacity, the sand will remain dangerously hot for a long time after the power is cut.
- Cleanliness: Ensure the sand is dry and free of organic contaminants, which could smoke or ignite at these temperatures

- Clean, dry silica or quartz sand is recommended. (Sand has to be arranged by the user)
- Bottom vs. Surface: Sand is an excellent insulator, which creates a significant temperature difference between the bottom and the surface
- The sand at the very bottom of the bath may reach 300 °C while the surface remains significantly cooler .
- Vessel Placement: To achieve 300 °C heating, the reaction vessel must be buried deep into the sand, near the heating source

Technical Specification

- Sand bath chamber size : 30 x 20 x 10 cm (L x W x H)
- Chamber volume : 06 Liters
- Outer Chamber : Fabricated from mild steel, which is powder coated in attractive shades
- Top heating plate : Made up of stainless steel Material
- Temperature Range : 50 °C to 300 °C (Higher range on request)
- Temperature Control : Microprocessor based Auto Tune PID Digital temp controller cum indicator
- Temperature sensor : RTD (Pt 100) sensor
- Temperature Display : LED Display
- Temperature Accuracy : ± 1 °C
- Control panel : Equipped with mains switch, control switch, pilot lamp and Digital Controller
- Heating : Heated by open type heaters placed at the bottom of the inner chamber
- Operation on : 230 Volts, 50 Hz, Single Phase, AC Supply
- Overall Dimensions : 330 X 430 X 365 (L X W X D) mm
- Weight (Approx.) : 15.500 Kg

Corporate Office:

A-1301, BVR Ek, Opp. Hotel Inder Residency,
Near Gujarat College, Ellisbridge,
Ahmedabad-380 006.
Phone : (079) 6621 1234
Website : www.eieinstruments.com

Calibration Laboratory:

Survey No. 1098, Nr. Mahakali Temple,
Opp. Primary School, Village Kubadthal,
Tal-Daskroi, Ahmedabad-382 433.
Phone: (079)-3520 8300
E-mail : info@eieinstruments.com

GST TIN No. : 24AABCE4018L1ZU
I.T. PAN No. : AABCE 4018L
IEC No. : 0813026598 Dt. 04/03/2014
CIN No. : U29199GJ2004PTC045078
TAN No. : AHME00572F